




BiTrac™ NIV

Non-invasive Ventilation masks



Critical Care ▪ Pulmonary®



BiTrac™ NIV

The BiTrac NIV product family features a variety of interfaces that are designed to increase patient comfort and compliance in respiratory distress scenarios. The full face and nasal masks are highly adjustable, and their seal's design minimise pressure exerted on the bridge of the nose.

Our range of interchangeable elbows greatly reduces waste, as the mask's elbow can be easily removed and replaced when transitioning patients between standard ventilators and bilevel machines.



Interchangeable elbows

Code	Description	Box Qty.
313-9599	BiTrac Select interchangeable standard non-vented elbow	10
313-9598	BiTrac Select interchangeable non-vented elbow with anti-asphyxiation valve	10
313-9593	BiTrac Select interchangeable vented elbow with anti-asphyxiation valve	10
313-9595	BiTrac Select interchangeable non-vented elbow with nebuliser port	10
313-9594	BiTrac Select interchangeable non-vented elbow with anti-asphyxiation valve and nebuliser port	10
313-9597	BiTrac Select interchangeable non-vented elbow with bronchoscopy port	10
313-9596	BiTrac Select interchangeable non-vented elbow with anti-asphyxiation valve and nebuliser port	10



BiTrac Select™ NIV full face mask with interchangeable standard non-vented elbow

Code	Description	Box Qty.
313-9509	BiTrac Select NIV full face mask with interchangeable standard non-vented elbow, extra large adult	10
313-9508	BiTrac Select NIV full face mask with interchangeable standard non-vented elbow, large adult	10
313-9507	BiTrac Select NIV full face mask with interchangeable standard non-vented elbow, medium adult	10
313-9506	BiTrac Select NIV full face mask with interchangeable standard non-vented elbow, small adult	10



BiTrac Select NIV full face mask with interchangeable non-vented elbow with anti-asphyxiation valve

Code	Description	Box Qty.
313-9504	BiTrac Select NIV full face mask with interchangeable non-vented elbow with anti-asphyxiation valve, extra large adult	10
313-9503	BiTrac Select NIV full face mask with interchangeable non-vented elbow with anti-asphyxiation valve, large adult	10
313-9502	BiTrac Select NIV full face mask with interchangeable non-vented elbow with anti-asphyxiation valve, medium adult	10
313-9501	BiTrac Select NIV full face mask with interchangeable non-vented elbow with anti-asphyxiation valve, small adult	10



BiTrac Select NIV full face mask with interchangeable vented elbow with anti-asphyxiation valve

Code	Description	Box Qty.
313-9514	BiTrac Select NIV full face mask with interchangeable vented elbow with anti-asphyxiation valve, extra large adult	10
313-9513	BiTrac Select NIV full face mask with interchangeable vented elbow with anti-asphyxiation valve, large adult	10
313-9512	BiTrac Select NIV full face mask with interchangeable vented elbow with anti-asphyxiation valve, medium adult	10
313-9511	BiTrac Select NIV full face mask with interchangeable vented elbow with anti-asphyxiation valve, small adult	10

**BiTrac Select™ NIV full face mask with interchangeable non-vented elbow with anti-asphyxiation valve and nebuliser port**

Code	Description	Box Qty.
313-9518	BiTrac Select NIV full face mask with interchangeable non-vented elbow with anti-asphyxiation valve and nebuliser port, extra large adult	10
313-9517	BiTrac Select NIV full face mask with interchangeable non-vented elbow with anti-asphyxiation valve and nebuliser port, large adult	10
313-9516	BiTrac Select NIV full face mask with interchangeable non-vented elbow with anti-asphyxiation valve and nebuliser port, medium adult	10
313-9515	BiTrac Select NIV full face mask with interchangeable non-vented elbow with anti-asphyxiation valve and nebuliser port, small adult	10

**BiTrac™ NIV nasal mask with standard non-vented elbow**

Code	Description	Box Qty.
313-9016W	BiTrac NIV nasal mask with standard non-vented elbow, extra large adult	10
313-9015W	BiTrac NIV nasal mask with standard non-vented elbow, large adult	10
313-9014W	BiTrac NIV nasal mask with standard non-vented elbow, medium adult	10
313-9013W	BiTrac NIV nasal mask with standard non-vented elbow, small adult	10

Replacement headgear

Code	Description	Box Qty.
313-9026	Replacement headgear for BiTrac full face and nasal masks	10
313-9027	Replacement CAP headgear for BiTrac full face and nasal masks	10

BiTrac, BiTrac Select and Omni-Clip are manufactured by Pulmonary Inc, USA and distributed by Intersurgical Ltd

Please think before you print Save energy and paper. If you must print this information sheet please print it double sided.

IS1.62 ▪ Issue 4 09.23



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com

UK ▪ Australia ▪ Benelux ▪ Brazil ▪ Canada ▪ China ▪ Colombia ▪ Czech Republic ▪ Denmark ▪ France ▪ Germany ▪ India ▪ Ireland ▪ Italy
Japan ▪ Lithuania ▪ Malaysia ▪ Mexico ▪ Norway ▪ Philippines ▪ Portugal ▪ Russia ▪ South Africa ▪ Spain ▪ Sweden ▪ Taiwan ▪ Turkey ▪ USA